

Astronomy to medical imaging: Blackford Analysis

Alan Heavens

Institute for Astronomy, University of Edinburgh

with Ben Panter, Rob Tweedie, Mark Bastin, Will Hossack, Keith McKellar, Trevor Whittley



ROE Workshop 13th October 2010



Blackford Analysis

- Blackford Analysis Ltd. is a spinout company from the University of Edinburgh that builds processing solutions using MOPED and other technologies
- MOPED is a patented technology which speeds up computation involving large datasets




Astrophysics Data Deluge

- ✦ Astronomy Datasets
- ✦ 1.4bn pixels every 30 sec
- ✦ > 1 TB / night image data
- ✦ Multi-object spectrographs
- ✦ Analysis Requires:
 - ✦ Accurate answers to fitting of complex models



Examples

- Galaxy spectra: flux measurements are sum total of starlight from stars of given **age** (simplest)






- 2 parameters = age, mass

Inverse problems

- ✦ **Parameter Estimation:** given some data, and a model, what are the most likely values of the parameters, and what are there errors?
- ✦ Astronomical problems often have
 - ✦ **LARGE DATASETS** and
 - ✦ (fairly) **LARGE PARAMETER SPACES**


Example: cosmology

- **Model:** Big Bang theory
- **Parameters:** Expansion rate, density of ordinary matter, density of dark matter, dark energy content... (around 15 parameters)



Dealing with large parameter spaces

- ✦ Don't explore it all - generate a chain of random points in parameter space
- ✦ Most common technique is MCMC (Markov Chain Monte Carlo)



Variants: Hamiltonian
Monte Carlo, Nested
Sampling, Gibbs
Sampling...

Large data sets


- ✦ What scope is there to reduce the size of the dataset?
- ✦ Why? Faster analysis
- ✦ Can we do this without losing accuracy?
- ✦ Depending on where the information is coming from, often the answer is yes.

MOPED weighting vectors

- ✦ MOPED automatically weights each data point in an optimum way, creating a **vastly** smaller dataset.
- ✦ In many cases, *the errors from the compressed dataset are no larger than those from the entire dataset*
- ✦ It is **NOT** obvious that this is possible


How does MOPED work?

Carefully constructed compression step allows rapid exploration of problem.



12D registration example: **1,000,000** combinations / second

Planets Around Distant Stars



Brute force: 3 years CPU
MOPED: 1 day CPU (1000x speed up)




10^{10}
combinations

light



Time mod period

Medical imaging: registration





Stroke lesion


Image Distortions

- ✦ MRI scans:
512x512x100 voxels
2.6 x 10⁷ data


- ✦ Affine distortions: **12 parameters**



Registration in MR / CT (implemented)





Registration takes **1 second**, down from **10 minutes** without **MOPED**




Transform:

[x, y, z] translations




[xy, yz, zx] rotations

[x, y, z] shears



[x, y, z] scales



2 x [512 x 512 x 100] volumes

12D => massive number of possible solutions

52 Million Voxels


What kind of problem is best?

Data

10010111010101000101111010100001010101000101001
 00101110101000000101000101000101000101110101000
 0101101010000001010001010001010001011101010001
 011110101000010101010001010001011101010001
 01000101010001011101010001011101010001
 01000101000

10010111010101000101111010100001010101000101001
 00101110101000000101000101000101000101110101000
 0101101010000001010001010001010001011101010001
 011110101000010101010001010001011101010001
 01000101010001011101010001011101010001
 01000101000

MOPED




Model

10010111010101000101111010100001010101000101001
 00101110101000000101000101000101000101110101000
 0101101010000001010001010001010001011101010001
 011110101000010101010001010001011101010001
 01000101010001011101010001011101010001
 01000101000


10010111010101000101111010100001010101000101001
 00101110101000000101000101000101000101110101000
 0101101010000001010001010001010001011101010001
 011110101000010101010001010001011101010001
 01000101010001011101010001011101010001
 01000101000

Large dataset
 Currently too slow
 Fitting a model


Planets Around Distant Stars



Linear Registration in MR / CT




Where could MOPED be applied?





Faster
Processing?

Smaller
Hardware?

Bigger
Datasets?



Hyperspectral Imaging



Acquire 512x512 spectra from hyper-spectral camera, each with 512 wavelengths. (Data alone is **134Mpix**)

Use MOPED **to fit spectral features on camera**, returning parametric fit for each pixel, vastly reducing data size and allowing further analysis

Problem: Spectral lines alone cannot distinguish

Answer: Fit whole spectrum?

-> Distinguish between chemical explosives and fuel

Determine spectral elements of interest – foliage, Semtex, chemical explosive, fuel, earth, rock & brick – MOPED filters to combine slices to reveal presence of spectra. Potentially eliminate real world lighting variations (reflections etc.) by incorporating as nuisance parameters

Original Datacube: **134Mpix**

Processed Data: **0.26Mpix** / parameter

10 'filters' – **50x** fewer data to transmit

Processing Speedup: **512/N_{parameters}**

Security Scanning



X-ray / THz scan identifies 30 objects in luggage

MOPED used to immediately compare each object to library of known threats in every orientation

Go/no-go advice given, or alerts operator for follow up physical examination

Multi-frequency scans allow spectral analysis of liquids
MOPED enables fast spectral fitting (see prev. slide)

Summary



- Blackford Analysis is using solutions designed for dealing with large Astronomical Datasets for rapid and accurate analysis of data
- Applications in medical imaging, security etc. - anywhere where speed and accuracy are important
- Arose from blue-skies research

Contact: Dr. Ben Panter

ben@blackfordanalysis.com

www.BlackfordAnalysis.com